
Metro PowerPod™ System

DESCRIPTION

Tripp-Lite AC Power Supply for use with Metro iMed, Flexline and Starsys carts.

Input: 100 VAC to 240 VAC, 50 Hz / 60 Hz.

Output: 120 VAC / 60 Hz, 150 VA or 150 W.

TABLE OF CONTENTS

COPYRIGHT NOTICE.....	2
SAFETY INFORMATION.....	2
INSTRUCTIONS FOR USE.....	4
DESCRIPTION OF MAJOR COMPONENTS.....	7
METRO POWERPOD™ SYSTEM.....	7
DISASSEMBLY OF POWER SUPPLY FROM CART.....	10
REPLACEMENT PARTS LIST.....	11
WIRING DIAGRAM.....	12

COPYRIGHT NOTICE

This manual is copyrighted. All rights reserved. This manual may be printed for personal use only. This manual, whole or in part, may not be copied, photocopied, reproduced, translated, or reduced to any electronic medium or machine-readable form for distribution. This manual, whole or in part, may not be modified without prior consent, in writing, from InterMetro Industries Corporation.

**Copyright © 2019 by InterMetro Industries Corporation, 651 North Washington Street,
Wilkes-Barre, PA 18705, U.S.A.
Phone: 1-800-992-1776
<http://www.metro.com/support>**

InterMetro Customer Service: For all customer service related issues, or if you need technical assistance, please call our customer service department at:

1-800-992-1776	Americas
+31 76 58 77550	Europe
+9714 811 8286	Middle East
+65 6567 8003	Asia-Pacific

SAFETY INFORMATION

This section contains important safety and operating instructions for PowerPod™ System. Please read all instructions on the cart before putting into service.

- CAUTION:** User maintenance consists solely of cleaning and minor external adjustments. For safety purposes, all servicing must be performed by qualified service personnel only. For all service items, please contact InterMetro Customer Service.
- CAUTION:** Do not operate the PowerPod™ System if it has received a severe impact or has been otherwise physically damaged. Please have a qualified service person inspect PowerPod™ System for any performance or safety hazard prior to putting back into service.
- CAUTION:** To prevent damage to any cords or connectors when disconnecting, always grasp and pull by the connector and not the cord. Do not operate the PowerPod™ System with damaged cords or connectors. Please replace any damaged component immediately. Contact InterMetro Customer Service for replacement parts and service.
- DANGER:** PowerPod™ Systems are not for use in hazardous (classified) locations. Do not use or recharge the power supply battery in oxygen enriched areas, areas where flammable anesthetics are used or stored, or any other hazardous, classified location.

NOTICE THE RECEIVER OF THIS PRODUCT IS RESPONSIBLE FOR FREIGHT DAMAGE CLAIMS.

- This equipment must be opened immediately for inspection.
- All visible damages must be reported to freight company within **48 hours**, and must be noted on freight bill at time of delivery.

- Concealed damage is your responsibility - you must advise carrier of any loss or damage within **15 days** after receipt of material.
- If there is damage, retain original packaging for inspection.

FCC Information:

The power supply (Tripp Lite HC150ATD) complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

Regulatory Compliance

The power supply (Tripp Lite HC150ATD) was tested to, and complies with, the following:

ANSI/AMMI ES60601-1
CAN/CSA-22.2 No. 60601-1
IEC 60601-1-2
IEC/EN 60601-1-2
Annex 1-1 (KN 60601-2)
FCC 15B
CAN ICES-3 (A)/NMB-3(A)
UL 1778
CSA C22.2 No. 107.3-14
NFPA 99

INSTRUCTIONS FOR USE

The Metro PowerPod™ System is designed to provide AC voltage output power from either utility provided AC input voltage or from inverting DC power from a battery. It may be configured for either a Lithium Iron Phosphate (Li-Ion) battery or a sealed lead acid (SLA) battery (SLA availability varies by cart product line.) To ensure proper operation, the following guidelines should be observed:

General:

1. The power module can support continuous AC loads up to a maximum of 150 VA or 150 W.
2. The power module input is rated 100 VAC to 240 VAC, 5.4 - 2.1 A, 50 Hz / 60 Hz (operating range from 90 VAC to 265 VAC and from 45 Hz to 66 Hz).
3. The power module output voltage is fixed at 120 VAC +/-5% at 60 Hz +/-0.3 Hz.
4. The power module is connected to the remote user interface module using a straight-through, six wire telephone type cable equipped with RJ-12 style connectors. The remote user interface module includes an audible alarm to indicate a low state of charge of the battery. The loudness of the audible alarm may be adjusted to loud, soft, or silent by setting the Alarm DIP switches on the power module. The power module is shipped with the switches set to soft mode.
5. The power module has internal over temperature protection. If the internal temperature exceeds the over temperature protection threshold, the power module will shut down (no VAC output and no charging). If the temperature drops below the threshold, charging will resume but the VAC output will remain off until the ON switch on the remote user interface module is pressed.
6. The power module has a charging fold back feature to prevent over temperature conditions. If the internal temperature of the power module becomes higher than normal while charging at 20 A, it will fold back the maximum charging current to 9.9 A, until such time as the power module cools to normal operating temperature.
7. When using a Li-Ion battery, a communication cable is required from the power module BATT connector to the battery communication connector. This allows the battery to communicate its operating parameters to the power module, including state of charge and over temperature conditions (the power module will interrupt current flow into or out of a battery if a battery reports an internal temperature alarm).
8. When using a SLA battery, a communication cable is not required. An optional temperature sensor is provided and is connected from the BATT connector to the negative terminal of the battery to provide added fault protection. The sensor will allow the power module to interrupt current into or out of a battery if a temperature is sensed that exceeds the over temperature threshold or if the sensor becomes disconnected. The sensor must be enabled by setting the Thermistor DIP switch to ON (set at factory).

Battery:

1. When changing battery type, power off the power module first, and then the battery selector DIP switch can be set to SLA or LI-ION.
2. A Li-Ion battery provides its % capacity level to the power module through the battery communication cable. The power module will calculate the % capacity of a SLA battery by measuring battery voltage and current.
3. The power module includes a low voltage cutoff (LVC) function that will remove any load on the battery when the battery is depleted to zero % reported capacity (Li-Ion) or 10.5 V (SLA).
4. Following the initial battery installation or a battery replacement, the power module should be allowed to charge the battery to full capacity before starting a discharge cycle of the battery. When a Li-Ion battery is first connected to a power module, there may be a delay of up to 60 seconds before full charging current is delivered to the battery.
5. If a Li-Ion battery has a voltage less than 10.8 V – 11.5 V, then the charger will start in recovery mode with a 1 A charging current. If the battery voltage recovers above the threshold within a timeout period, then full charging current will be resumed. If not, the battery will be considered defective and the charger will shut down, but 120 VAC output will still be provided.
6. If a SLA battery voltage has a voltage less than 9.6 V, then the charger will start in recovery mode with a 1 A charging current. If the battery voltage recovers above the threshold in less than 120 minutes, then full charging current will be resumed. If not, the battery will be considered defective and the charger will shut down, but 120 VAC output will still be provided.

Installation:

1. **Caution:** To assure optimum performance and battery life, set the battery selector DIP switch to the proper battery type before connecting the battery.
2. **Caution:** To prevent over temperature issues for the battery, the wiring, or the power module, the DC power cable to the

battery must be tightly secured to the battery terminals.

3. **Caution:** To prevent possible damage to the power module when using a Li-Ion battery, make the battery DC power connections before making the battery communication connections (and reverse the procedure when disconnecting the battery).
4. **Caution:** Do not drill into or attempt to open any part of the power module housing, there are no user-serviceable parts inside. Do not attempt to use the power module if any part of it becomes damaged.
5. **Caution:** To avoid the risk of electrical shock, only connect this equipment to a supply mains that provides a protective earth connection.
6. **Caution:** The EMISSIONS characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.
7. **Caution:** Output voltage is only supported for a limited time, depending on battery size and load. Suitability of runtime on battery needs to be evaluated in the end use application.
8. **Caution:** All SIP/SOP (SIGNAL INPUT/OUTPUT PART) Connections are for SELV (Safety Extra-Low Voltage) Only.

Statement of Intended Use:

Caution: Not intended to power life support equipment.

Caution: Do not use this equipment within oxygen enriched atmospheres, or within 0.3 meters of a point at which oxygen enriched atmosphere is intentionally vented.

Environmental conditions:

Storage temperature: -15 to 50 Degrees C. (5 to 122 Degrees F.)

Operating temperature: 0 to 40 Degrees C. (32 TO 104 Degrees F.)

Humidity: 5% to 95% RH (non-condensing)

Symbols and Markings:

Class I = Protective Earth required.

= Refer to Instruction (Application Guide)

= Alternating Current.

IPXO = No Ingress Protection.

NOTICE For batteries recycling/disposal in the United States and Canada, call2Recycle.

NOTICE For other countries, batteries must be recycled / disposed of in compliance with local regulations.


NOTICE

Contact InterMetro service for product repair or battery replacement. Use of unauthorized service or batteries will void warranty, damage equipment, and / or affect performance. C05-1100

METRO POWERPOD™ SYSTEM — INSTRUCTIONS FOR USE

DANGER: Metro PowerPod™ Systems are not for use in hazardous (classified) locations. Do not use or recharge the power supply battery in oxygen enriched areas, areas where flammable anesthetics are used or stored, or any other hazardous, classified location.

Shipment of Power Supply

Depending on the battery chemistry, the preparation and marking of the shipment are different.

Lithium Iron Phosphate Battery (Li-Fe)

The shipping container must have the following marking:

Class 9 Shipment:

Outside the USA in any shipping mode and air or sea transport inside the US, the power supply must be shipped as fully regulated Class 9 (Hazardous Material / Dangerous Goods).

In this shipping mode the batteries can be connected or disconnected. The shipper must be certified for Li-Ion Class 9 shipments.

The power supply must be shipped per UN3481 Lithium ion batteries contained in equipment on cargo aircraft only.

SLA

The SLA battery is not considered a Hazardous Material (Dangerous Goods) and there is no restriction on shipping the SLA batteries.

DESCRIPTION OF MAJOR COMPONENTS

Depending on your order, your PowerPod™ System will be delivered in a specific configuration. The configuration may employ the following major components:


- Power Supply
- Remote User Interface (RUI)
- Power (outlet) strip
- Power Cord

Power Supply

- This Power Supply is a fully automatic power supply charger system with an AC output voltage of 120 Volts, 60 Hz
- When plugged into an AC outlet, it supplies power to the equipment while also charging the battery.
- When unplugged from the outlet, the power supply switches automatically to supply the equipment from the internal battery.
- The power supply is designed for continuous operation and service is not interrupted by plugging in or unplugging the power supply.

Remote User Interface (RUI)

- To monitor the Mode status and charge level of the power supply battery. See Figure.
- The RUI is located at the back center of cart and gives the user the current battery status at eye level.

- The RUI is connected to the power supply by a cable from the power supply.

Power (outlet) strip

- Electrical equipment on the cart (computer, monitor, etc.) is powered by plugging them into the power strip located on the rear of the cart. **NOTE: Total device power requirements must not exceed 150 W.** Outlet Covers are provided with North American units and are to be installed in unused outlets.

Power Cord

- The Power Cord is selected on the basis of Domestic or International Application. For International, the AC In Power Cord is selected by Configurator based on the Country Specified

METRO POWERPOD™ SYSTEM

Connect Power Supply

The unit is shipped with the battery cable unplugged from the power supply. The cable must be plugged in and the battery charge level checked before using the unit.

1. Ensure unit is unplugged from wall outlet.
(If the unit is ON, turn OFF by pushing and holding the “MODE” button for 2 seconds until the display goes blank.)

METRO POWERPOD™ SYSTEM — INSTRUCTIONS FOR USE

- Unplug the power strip cord and the display communication cable from the PowerPod™ drawer at the rear of the cart (See **Image 1**).


Image 1: Rear of PowerPod™

- Loosen the lock screw and remove any security seal (if present) from the front of the PowerPod™ drawer (See **Image 2**).


Image 2: Front of PowerPod™

- Slide the drawer latch to the left to unlock the PowerPod™ drawer.
- Open drawer and connect the battery / power supply cable (See **Image 3**). (Match red to red and black to black)


Image 3: Battery/Power Supply Cable

- Close drawer slide, latch to the right and tighten screw and replace security seal (if desired).
- Reconnect outlet strip power cable and communication cable to the PowerPod™ at the rear of cart (See **Image 1**).

Charge Battery

Battery must be charged for a minimum of 4 hours prior to initial startup.

- Plug the PowerPod™ power cord directly to a wall outlet to charge battery.
- After the PowerPod™ reaches full charge, turn the power system ON by pressing and holding the Mode button on the fuel gauge until the display comes on (See **Image 4**).


Image 4: Full Charge

- Devices may now be plugged into the outlet strip located on the rear of the cart.
To maintain Medical Grade requirements, install provided outlet covers in unused outlets (North American domestic carts only).

Note: It is recommended that the cart be plugged in and charging when not in use. The power system is equipped with a smart power supply that will prevent the battery from being overcharged when plugged in for extended periods of time.

Beeper Volume

The beeper has three volume settings: LOUD, SOFT, or OFF. The default setting is SOFT.

The beeper volume control is located on the power supply in the drawer.

- Loud = DIP switches 2 & 3 in UP position
- Soft = Dip Switch 2 in UP position, DIP Switch 3 in DOWN position
- Off = DIP switches 2 & 3 in DOWN position

Electrical Input: 100 VAC to 240 VAC, 5.4 - 2.1 A, 50 Hz / 60 Hz

Electrical Output: Loads up to a maximum of 150 VA or 150 Watts; 120 VAC, 60 Hz

Note: Total device power requirements must not exceed 150 W.

When first plugged in the system will be in Standby Mode and not be sending power to the outlet strip. To enable power, press the “MODE” button for 2 seconds until the runtime minutes changes to 999.

The power system has three operating modes:

1. Standby Mode (When plugged into a live AC outlet)

This is the mode when a unit is first plugged into the wall receptacle. Denoted by **flashing** AC input icon  and 0 minutes runtime shown.

- Unit is charging the battery only; AC output is OFF
- There is no power being sent to the outlet strip.

To turn system off, unplug from AC outlet.

To **enable power** (Charge Mode), press and hold Mode button for 2 seconds.

2. Charge Mode (When plugged into a live AC outlet)

This is the mode that provides power and is activated from Standby Mode by holding the Mode button for 2 seconds.

Denoted by a **solid** AC input icon  and 999 minutes runtime shown.

- Unit is charging the battery and AC Output is ON
- Power is now being sent to the outlet strip.

To **disable power** (Standby Mode), press and hold Mode button for 2 seconds.



To turn system off, disable power then unplug from AC outlet.

3. Battery Mode (When cart is **not** plugged into a live AC outlet)


This mode is switched to automatically from Charge Mode, when the cart is unplugged from a live AC outlet.

Denoted by the **absence** of the AC input icon and the estimated runtime in minutes is shown (i.e. will be a number other than 0 or 999).


- This occurs when the unit is ON, but not plugged into a [live] AC outlet.
- The battery is not being charged.
- When battery capacity is less than 10%, the alarm will beep (if sound set to soft or


MODE Button

Standby Mode Indication

Charge Mode Indication

Battery Mode Indication

METRO POWERPOD™ SYSTEM — INSTRUCTIONS FOR USE

loud) and the display will flash every 3 seconds. Both can be stopped by pressing the “Mute” button.

- When battery capacity is 5% or less, the alarm beeps (if sound set to soft or loud) and the display blinks every 3 seconds. The audible alarm can be **muted for one minute** by pressing the “MUTE” button. The display will continue to blink. The battery should be charged at this time to prevent unwanted shutdown.
- To turn system off (**disable power**), press and hold Mode button until the display goes blank.
- To turn system on (**enable power**), press and hold Mode button until the display lights up.

DISASSEMBLY OF POWER SUPPLY FROM CART

WARNING: Before disassembly, power cord must be disconnected from the wall outlet.

NOTE: Disconnect all wire harness and electrical cords from power supply.

- Loosen the lock screw and remove any security seal (if present) from the front of the PowerPod™ drawer (See **Image 2**).


Image 2: Front of PowerPod™

- Slide the drawer latch to the left to unlock the PowerPod™ drawer.
- Open drawer and disconnect the battery / power supply cable (See **Image 3**).


Image 3: Battery/Power Supply Cable


REPLACEMENT PARTS LIST

Certain replacement parts are based on the electrical requirements of the locality.
Please confirm desired plug type or other considerations when ordering.


Index #	RP Part Number	Description
1	RP-PPSLIDES	RP TRAY SLIDES, POWERPOD
2	RP-PPTRAYASSY	RP TRAY ASSY, POWERPOD
3	RP-PPNEMAINPCORD	RP NEMA/TYP E B INPUT CORD, POWERPOD
4	RP-PPNEMAOUTCORD	RP NEMA/TYP E B OUTPUT CORD, POWERPOD
5	RP-PPXPINPCORDEF	RP EXPORT TYPE E/F INPUT CORD, POWERPOD
6	RP-PPXPINPCORDG	RP EXPORT TYPE G INPUT CORD, POWERPOD
7	RP-PPXPINPCORDI	RP EXPORT TYPE I INPUT CORD, POWERPOD
8	RP-PPXPINPCORDL	RP EXPORT TYPE L INPUT CORD, POWERPOD
9	RP-PPXPOUTCORD	RP EXPORT OUTPUT CORD, POWERPOD
10	RP-PP30AFUSES	RP 30 A FUSES (5/PKG), POWERPOD
11	RP-PPLIFEBATT	RP LIFE BATTERY, POWERPOD
12	RP-PPSLABATT	RP SLA BATTERY, POWERPOD
13	RP-PPPWRSUPPLY	RP POWER SUPPLY, POWERPOD
14	RP-PPPCUKIT	RP PCU (RUI) KIT, POWERPOD
15	RP-PPPCUCABLE	RP PCU (RUI) CABLE, POWERPODD
16	RP-PP4PWRSTRIP	RP 4 OUTLET POWERSTRIP, POWERPOD
17	RP-PPXPPWRSTRIP	RP 6 OUTLET EXPORT POWERSTRIP, POWERPOD
18	RP-PPBATTCABLE	RP MAIN BATTERY CABLE, POWERPOD
19	RP-PPDRAGCABLE	RP DRAG CABLE, GROUND, POWERPOD
20	RP-PPLATCH	RP LATCH WITH HARDWARE, POWERPOD
21	RP-PPMNTRAIL48	RP LMED400/800 MOUNTING RAILS W/HW, POWERPOD
22	RP-PPMNTRAIL6	RP LMED600 MOUNTING RAILS W/HW, POWERPOD
23	RP-PPMNTRAILFLEX	RP FLEXLINE MOUNTING RAILS W/HW, POWERPOD
24	RP-PPMNTRAILSTAR	RP STARSYS MOUNTING RAILS W/HW, POWERPOD
25	RP-PPLIFTSTOP	RP BATT LIFTER, STOP, W/HW, POWERPOD

Power cord replacement number based on type of plug.


WIRING SCHEMATIC FOR DOMESTIC WITH SLA BATTERY

WIRING SCHEMATIC FOR INTERNATIONAL WITH LiFePO4 BATTERY


WIRING SCHEMATIC FOR INTERNATIONAL WITH SLA BATTERY

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Copyright © 2019 par InterMetro Industries Corporation, 651 North Washington Street, Wilkes-Barre, PA 18705, États-Unis