

Advancing the way you care.



Flexline® Telehealth

The Flexline Telehealth Cart provides more room to work, with a generous worksurface, large drawers and the addition of PowerPod™ to power a computer, monitor and webcam. It is the perfect all in one solution to contain the needed patient care supplies plus a connection to skilled care, no matter where the physician might be.

MICROBAN  **24/7 Protection**

Long lasting power & protection when you need it the most.

Flexline carts are constructed with a polymer body and all polymer drawers that won't rust, dent, chip or peel. The polymer is molded with Microban® antimicrobial product protection built in. Microban® antimicrobial protection works 24/7 for an added level of protection against bacteria, mold and mildew and makes products easier to clean and keep cleaner in between cleanings.

METRO POWERPOD™

Computerized carts are equipped with Metro's PowerPod, featuring a universal power supply that accommodates a long running LiFe battery solution that provides a lower cost of ownership. PowerPod is an AC/ "Plug-& Play" solution designed to get the clinician and tools they need—where they are needed...Next to the patient.



Battery Chemistry	LiFe <i>Lithium Iron Phosphate</i>
<i>Battery Capacity</i>	576 Wh
<i>Available Power*</i>	490 Wh
<i>Run Time @ 35 Watts</i>	14 hours
<i>Battery Life (@ 60%)</i>	3,000+ cycles
Power Supply	<i>Universal</i>
<i>Recharge Time (max.)</i>	2.5 hours
<i>Charging Status Indicator</i>	Continuous readout on control unit display
<i>Input Voltage</i>	100-240 VAC, 5.4-2.1 A, 50 Hz / 60 Hz
<i>Output Voltage</i>	120 VAC, 60 Hz
<i>Maximum Power Draw</i>	150 W
<i>Regulatory</i>	UL 1778, CSA C22.2 No. 107.3-05; CAN/CSA-22.2 No. 60601-1-08; ANSI/AMMI ES60601-1:2005; Complies with NFPA 99; Marks: CE, FCC and EAC

* Power to device after converting DC battery power to AC and accounting for depth of discharge protection.

Visit www.metro.com or call **1-800-992-1776**.
 LO7-175 7/20
 © 2020 InterMetro Industries Corporation, Wilkes-Barre, PA 18705

