

# Which cart should I choose?

While there are many different types of procedures, from Arthroscopic to Vertebroplasty, they typically have one thing in common, the need to accommodate specialty equipment. The procedure defines the type of equipment and accessories required which is why a procedure cart should have a wide array of accessories that can handle supply and equipment needs. A good procedure cart will also allow for easy adaptability in the field. With features like these, plus an optional on-board power source and expandable workspaces, Metro's medical carts give you more flexibility and adaptability for your budget, giving you the power to do more.



GOOD



Basix Plus  
No. **MBP3210TRAUMA**

**Traditional design  
& construction.**



- Lighter weight hybrid construction
- Clean corner attachment
- Configure without penalty



BETTER



Starsys  
No. **SXRSARTH**

**Modern design &  
highly configurable.**



- Lightweight polymer construction
- Antimicrobial product protection
- Small footprint, expandable workspace
- Optional advanced wireless security



BEST



Flexline  
No. **FLWENDO**

**Modern design with  
expansive capacity.**



- Larger drawers and workspace
- Antimicrobial product protection
- Lightweight polymer construction
- Optional advanced wireless security



## CLEANER BY DESIGN™

- Microban® antimicrobial protection helps to keep products "cleaner between cleanings."
- Advanced polymer and other proprietary finishes provide corrosion resistance.
- Smooth, rounded corners to allow for easier cleaning.

\*Microban protection is not designed to protect users against disease causing microorganisms.