

Which cart should I choose?

Balancing storage capacity with available floor space and features versus costs are important considerations when choosing an isolation cart. However, one concern less thought about is the identification of an isolation room without undo alarm to patients, family members, and visitors. Yellow is the traditional color identified with isolation. It can be as bold as an entire cart or as subtle and direct as a drawer front.


GOOD



Basix Plus
No. **MBP0120ISO**

**Traditional design
& aesthetics.**

- Lighter weight hybrid construction
- Clean corner attachment
- Configure without penalty


BETTER



Flexline
No. **FLISO3**

**Modern design
with room to move.**

- Spacious storage and workspace
- Lightweight polymer construction
- Antimicrobial product protection
- Configure to preference


BEST



Starsys
No. **SXRSISO1**

**Modern design
with room to store.**

- Small footprint
- Lightweight polymer construction
- Antimicrobial product protection
- Configure to preference



CLEANER BY DESIGN™

- Microban® antimicrobial protection helps to keep products "cleaner between cleanings."
- Advanced polymer and other proprietary finishes provide corrosion resistance.
- Smooth, rounded corners to allow for easier cleaning.

*Microban protection is not designed to protect users against disease causing microorganisms.