

Which cart should I choose?

An anesthesia cart must be durable and secure to withstand the demands associated with daily use and ensure medications remain accessible only to authorized staff members. Good organization is key for easy identification and quick retrieval of medications, tools, and supplies. Additionally, there must ample storage to assure that all essential items have a defined location in or on the cart.



GOOD



Basix Plus
No. **MBP3210ANES1**

Traditional design that matches existing fleet. 

- Lightweight hybrid construction
- Clean corner attachment
- Configure without penalty



BETTER



Starsys
No. **SXRWSANES2**

Modern design with expandable workspace. 

- Advanced wireless security
- Lightweight polymer construction
- Antimicrobial product protection
- Expandable workspace



BEST



Flexline
No. **FLWANES1**

Modern design with advanced capacity. 

- Advanced wireless security
- Integrated narcotics drawer
- Integrated side storage
- Antimicrobial product protection
- Lightweight polymer construction



CLEANER BY DESIGN™

- Microban® antimicrobial protection helps to keep products "cleaner between cleanings."
- Advanced polymer and other proprietary finishes provide corrosion resistance.
- Smooth, rounded corners to allow for easier cleaning.

*Microban protection is not designed to protect users against disease causing microorganisms.